SERACONTROL N

REF 5010605R 6x5 ml 
 5011805R 18x5 ml
 For in vitro medical device
 Per uso diagnostico in vitro

Uso

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Use

SERACONTROL is a multiparametric control which intended use is control for clinical chemistry assay.

Summary

Lyophilized control serum prepared from human serum. The concentration and activities of the components have been adjusted for a wide range of automatic analyzers to ensure optimal control references of the clinical chemistry reagents.

Calibration values

Concentrations and activities of the components are lot-specific. The exact values are given in the data sheet enclosed in the kit.

Preparation

Open one vial of control avoiding the loss of lyophilizate.
 Reconstitute the lyophilized serum with exactly 5.0 ml of bi-distilled water.
 Close the vial and let it stand for 30 minutes out of bright light.
 Before use mix by inverting the vial.
 Do not shake the vial as the formation of foam should avoid.
 Ensure that all lyophilizate is dissolved.

Storage and stability

Prior to reconstitution, when stored at 2-8 °C, the control is stable until the expire date.

After the reconstitution the serum is stable:

8 hours	if stored at	+25° C
7 days	if stored at	+4° C
30 days	if stored at	-20° C

- The bilirubin stability, if stored protected from light exposure, is about 4 days at +2/8°C. Do not store at +15/25°C and do not freeze.
- The bicarbonate stability is 8 hours in the closed bottle and 1 hour in the open bottle.
- Alkaline phosphatase values in the reconstituted serum tend to rise during the stability period.
- For total and prostatic acid phosphatase the material should be stabilised by adding 1 drop (25-30 µl) of 0.7 M acetic acid solution to 1 ml of serum. After stabilisation total and prostatic acid phosphatase are stable for at least 2 hours at +25°C, 2 days at +4°C and 1 month when frozen once at -20°C.

Note

- Avoid direct light exposure, contamination and evaporation.
- In case of compliant or quality control request, refer to the lot number on the package or to the lot number of the single vial.
- Non properly reconstitution and / or storage can feed erratic results.

Precaution in use

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998).

 Human serum is obtained using only blood of donors tested by FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV.

However, no known test method can offer complete assurance that product derived from Human blood transmit infectious diseases, this product should be handled as potentially infectious biological material with care according to good laboratory practice.


Waste management

Please refer to the local legal requirements.

Bibliography

Department of Labor, Occupational Safety and Health Administration 29 CFR Part 1910.1030 Fed. Register July 1, 1998:6:267-280
 Council Directive 90/679/EEC. Official journal of the European directives N°L374 from Dec.31,1990:1-12

Symbols

	Biological risk
	CE Mark (98/79 CE regulation)
	in vitro medical device
	Batch Code
	Use by
	Storage temperature limits
	Read instruction for use
	Gesan Production srl

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